


Zyxin Antibody (R35413)

Catalog No.	Formulation	Size
R35413-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7791
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Zyxin antibody is available for research use only.



Western blot testing of peripheral blood mononucleocyte lysate with Zyxin antibody at 0.1ug/ml. Predicted molecular weight: ~61kDa, routinely observed at ~82kDa.

Description

Additional name(s) for this target protein: HED2, ESP2; ZYX

Application Notes

Optimal dilution of the Zyxin antibody should be determined by the researcher.

Immunogen

Amino acids QQREKPRVQEKQH were used as the immunogen for this Zyxin antibody.

Storage

Aliquot and store the Zyxin antibody at -20oC.

References (1)